
**INSTALLATION MANUAL  
MANUEL D'INSTALLATION  
MANUAL DE INSTALACIÓN  
MANUALE D'INSTALLAZIONE  
EINBAUANLEITUNG  
MANUAL DE INSTALAÇÃO**

**HIDRAULIC CHAIR LIFT  
L'ELEVATEUR HYDRAULIQUE  
ELEVADOR HIDRÁULICO  
ELEVATORE IDRAULICO  
HYDRAULISCHER BEHINDERTENLIFT  
ELEVADOR HIDRÁULICO**

## 1. PACKAGING

The hydraulic chair lift is packed in three boxes which contain the following items:

- Handrails
- Guide cylinder assembly
- Seat, steps, anchors and other remaining components

## 1. OPERATION

The hydraulic chair lift basically consists of a hydraulic cylinder which operates with water at a pressure of between 3 and 5 bar. At 5 bar, loads of up to 115 kg can be lifted and this maximum load will reduce at lower pressures. The supply pressure should never exceed 5 bar. Each lifting cycle will use approximately 7 litres of water.

To enter the pool, the chair should be in the top position and once the user is seated, the lowering control should be operated. The chair will then start to descend by the venting of the water from the cylinder and as the chair lowers it simultaneously rotates 100° to allow the user to swim directly into the pool.

To leave the pool, the above procedure is reversed.

## 3. INSTALLATION AND ASSEMBLY INSTRUCTIONS

### 3a. ASSEMBLING THE LADDER

The position of the chair lift will be dependent on the deep of the pool (this installation is recommended in a pool with a minimum depth of 1,6m).


Fit the treads (29) between the two hand rails (23 & 30) as shown in the drawing, pass the bolts (31) through the holes in the hand rails and the threads, fit the washers (18) and tighten them down with the corresponding nuts (19).

Press on the footpads (22) to the bottoms of the handrails (23 & 30) to prevent the ladder from damaging the pool walls.

Assemble the adjustable supports as follows:

Completely screw the upper support tube (28) into the lower support flange (34), screw this assembly into lower support tube (27) and then fit the non threaded end of the entire assembly into the adjustable foot pad (26). This assembly should now fixed to the handrail as near as possible to the edge of the pool by using the upper support flange (35) in conjunction with the nuts (33), bolts (19) and washers (18) provided.


A special adjustable support kit can be requested for pools with coping stones or other special pools.

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
### 3b. ANCHOR INSTALLATION

Fit the two anchors to the handrails and fix them using the corresponding screws and washers.

Connect the earth cable to the anchors.

With the ladder assembled and correctly adjusted for level using the adjustable supports, embed the anchors into the concrete.

### 3c. ASSEMBLY OF THE LIFTING MECHANISM

Once the anchors have been installed, slide out the ladder from them and remove the ladder from the pool. Screw the valve support (44) to the seat using the relevant nuts, bolts and washers (41, 43 & 46). Connect the 2 elbows (5) to the valve (45) – Remembering to use PTFE tape on the threads – and then mount the valve onto the support (44) by using the thread nipple (42 – Again using PTFE tape on the threads). Connect the inlet hose (7) (0,9m) to the upper elbow and the outlet hose (47) (0,7m) to the lower elbow and retain the hoses in place by using jubilee clips (6).

Place the seat over the cylinder, after having first removed 3 of 4 screws which retain the upper guide plate (1).

Re-fit these screws and then tighten down the seat clamps using the various nuts, bolts and washers (18, 19, 48). Do not completely tighten down the clamps, as we want a bit movement in order to be able to orientate the seat.

Pass the bolts (25) through the handrail, concave seats (15) and spacers (17) and position the cylinder/seat assembly so that the guide rail (14) is towards the bottom of the ladder. Fix the cylinder in place by tightening down the nuts and washers (11, 12).

Fix the elbow (5) to the lifting cylinder (8), connect the hose to the elbow and tighten it down using the jubilee clip.



Connect the water inlet to the control valve (45), but do not allow water to reach the valve.

Remove the clip fitted to the lower section of the cylinder (this clip was only used to retain the cylinder during shipping and assembly) and lower the cylinder to its lowest position. Make sure that the seat does not touch the ladder.

Refit the ladder/hydraulic lift assembly in the anchors, ensuring that it is correctly level before fixing it in place.

Allow the water to reach the control valve (45) and open the valve so that the seat rises. Operate the control valve several times to check the correct operation of the lifting mechanism.

Position the seat as required and tighten down the clamp bolts (48) so that the seat remains firmly fixed in position.

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